SMARTeZ[®]

Exclusively Distributed by:

PROGRESSIVE

Elastomeric Infusion Pump

The SMART Solution for Chemotherapy Ambulatory Infusion

Enables Patients to Participate in Normal Daily Activities and Improve Quality of Life



MEDICAL



SMARTeZ® Elastomeric Infusion Pump

SMARTeZ® Chemotherapy-rated pumps are approved for use with chemotherapy medications. 5-Fluorouracil is commonly infused over a specific infusion duration. Please refer to the table on page 2 for example infusion durations that can be accurately achieved with over or under filling SMARTeZ® pumps.



SMARTeZ® Features:

- Components are not manufactured with natural rubber latex and are DEHP free
- Individually calibrated for accuracy
- Color-coded by flow rate to assist in filling accuracy and reduce the risk of medication errors
- Double cap to protect sterility during filling
- Ideal infusion pump for patients and pharmacy staff:
 - Durable, Lightweight & Portable
 - No batteries or programming required
 - Easy to use and simple to fill design for smoother and faster filling
- Automated manufacturing provides a consistent and reliable supply, diversified therapy options and customizable combinations
- Disposable and single use
- Chemotherapy rated pumps are designed for 5-FU infusions:
 - Designed with an amber-tinted outer covering
 - Flow restrictor element designed to reduce the risk of occlusion from drugs prone to precipitation

SMARTeZ® Pumps vs. Mechanical Pumps

Ambulatory pumps used to administer home infusion are no longer separately billable to Medicare.* SMARTeZ may offer cost savings when compared to use of mechanical pumps.

ZERO Mechanical Pump Cost of Ownership

ZERO Maintenance, Calibration, or Cleaning

ZERO Capital Equipment Inventory Management

ZERO Third Party Vendor Rental Payments

ZERO Nuisance Alarms

Mechanical Pumps Total Cost of Ownership

- ☑ Maintenance
- ✓ Inventory Mgmt
- ☑ Rental Payments
- ✓ Nuisance Alarms

A Pump for All Common 5-FU Infusions

Infusion Durations Commonly Dispensed with 5-Fluorouracil (SMARTeZ Over Fill and Under Fill Volumes to Achieve 5-FU Infusion Durations)									
Infusion Duration (hours)	Item #	Under/Over Fill Volume (ml)	Nominal Flow Rate (ml/hr)						
24	SE0004-120C	96	4						
24	SE0010-270C	240	10						
46	SE0002-100C	92	2						
46	SE0005-270C	230	5						
48	SE0002-100C	96	2						
48	SE0005-270C	240	5						
72	SE0005-270C	360	5						
96	SE0002-270C	192	2						
120	SE0002-270C	240	2						
144	SE0002-270C	288	2						
168	SE0002-270C	336	2						

[Pump Volume / Labeled rate of ml/hr = Infusion Duration Prescribed]

Example: 4000mg 5-FU over 46 hours | Pump SE0002-100C

[under filled to 92ml / 2ml/hr = 46 hour infusion]



The double cap design includes a disposable protective port cap. The winged design allows for easy removal and stands upright to protect the sterility of the color-coded fill port cap while filling the device. After filling the device reattach the color-coded port cap and discard the disposable winged cap.



^{*} Sourced from: http://www.cms.gov/Regulations-and-Guidance/Guidance/Transmittals/2016-Transmittals-Items/SE1609.html.

Ordering Information:



Туре	PMI Item #	Nominal Fill Volume (mL)	Nominal Flow Rate (mL)	Nominal Duration (Hrs)	Approx. Nominal Duration (Days)	Qty. per Case	Color Code
5-FU	SE0002-100C	100 mL	2 mL/h	50 Hrs	2 Days	12	Yellow
	SE0002-270C	270 mL	2 mL/h	135 Hrs	5.5 Days	12	Yellow
mo	SE0004-120C	120 mL	4 mL/h	30 Hrs	1 Day	12	Orange
Chel	SE0005-270C	270 mL	5 mL/h	54 Hrs	2 Days	12	Brown
	SE0010-270C	270 mL	10 mL/h	27 Hrs	1 Day	12	Light Green

Accessory: Elastomeric Carry Pouch



Item # PMIPOUCH Item # POUCH2XL

Visit us for more information and to learn more about our entire EPIC Medical product portfolio.



