

PMI Suture Guide

(10/12mm; 12mm-2 Holes; 12mm-4 Holes)



Instructions For Use



- ITEM NUMBERS:**
- PMISG10:** (10/12mm)
 - PMISG12B2:** (12mm – 2 Holes)
 - PMISG12B4:** (12mm – 4 Holes)

DESCRIPTION

Reusable, Non-Sterile PMI Suture Guides are packaged by style as singles:


- One (1) 10/12mm Guide
- One (1) 12mm Bariatric Guide – 2 Holes
- One (1) 12mm Bariatric Guide – 4 Holes

INDICATIONS FOR USE

The PMI Suture Guide is intended for use during endoscopic / laparoscopic surgery. It is to be used only by surgeons trained in endoscopic / laparoscopic surgery.

INSTRUCTIONS FOR USE

The PMI Suture Guide can be used in conjunction with the PMI Trocar Closure-Suture Grasper to close trocar wounds in the following manner (see figures below):



STEP 1

Insert PMI Suture Guide with holes aligned perpendicular to the incision. Push PMI Trocar Closure-Suture Grasper through PMI Suture Guide, fascia, muscle, peritoneum and into abdomen. Drop suture and remove PMI Trocar Closure-Suture Grasper.

STEP 2

Push PMI Trocar Closure-Suture Grasper through opposite side of PMI Suture Guide. Pick up suture.

STEP 3

Pull the suture up through the peritoneum, muscle, fascia & PMI Suture Guide.

STEP 4

Remove PMI Suture Guide. Tie suture to complete the closure.

PMI Suture Guide

(10/12mm, 12mm-2 Holes; 12mm-4 Holes)

ITEM NUMBERS: PMISG10 & PMISG12B2 & PMISG12B4

CARE AND HANDLING INSTRUCTIONS

INSTRUMENT CLEANING:

- A. Rinse with de-ionized water to remove any gross debris. Hand wash instrument using neutral pH detergent or enzymatic cleaner and a soft brush, being sure to flush all channels with the solution. Avoid the use of steel wool, wire brushes and/or abrasive detergents.
- B. Clean entire surface of instrument.
- C. Thoroughly rinse instrument with de-ionized water after cleaning to remove any residual debris or cleaning solution.
- D. Instruments that are to be steam sterilized should be lubricated with non-silicone antimicrobial, water soluble, lubricant per manufacturer's instructions. Do not rinse or towel dry. Do not use mineral oil, petroleum jelly or silicone spray, which may inhibit sterilization and cause build-up in the crevices of the instrument.

STERILIZATION:

Follow the hospital Infection Control Procedure. The instrument can be sterilized using a steam autoclave, heat sterilization or gas sterilization. Below are the recommended sterilization parameters. The use of "flash" sterilization is not recommended as it can shorten the life of the instrument.

Below are the recommended sterilization parameters:

Sterilizer	Exposure Temp.	Exposure Time	Minimum Dry Time
Pre-Vacuum (wrapped)	121° C (250° F)	20 min	20 min
	132° C (270° F)	4 min	20 min
	134° C (273° F)	3 min	15 min
Pre-Vacuum (unwrapped)	132° C (270° F)	4 min	
Gravity Steam (wrapped)	132° C (270° F)	18 min	

STORAGE:

Instruments must be completely dry and lubricated to prevent damage. Instruments should be stored in areas that provide protection from extremes in temperature and humidity.



Distributed by:
**Progressive Medical
Inc.**

997 Horan Drive
Fenton, MO 63026
314-961-5786

NON-STERILE

Non-Sterile

Rx

Professional Use Only



Read Instructions Before Use

Not Made With Natural Rubber Latex

CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician.