


Q2® MultiPort® from QUEST Medical, Inc.

Minimal complexity —  
**Maximum protection.**



MultiPort®  
9520



Typical (complex)

Q2 MultiPort (simplified)

### Six ports; one flow

Give your patients the quality care they deserve with Q2® split septum, needle-free MultiPort Manifold. In-line gravity check-valves provide one way flow without manipulation of stopcock handles. Quest's Q2 fluid delivery systems offer convenient designs, configured to simplify the delivery of medications.

### Designed for safety

- Straight, direct fluid pathways for optimum flushing
- Closed systems designed to reduce the chance of contamination, maintaining a sterile environment<sup>1</sup>
- Low residual volume/dead space
- Neutral fluid displacement when attaching and detaching
- Q2 check-valves ensure one-way fluid delivery of medication to the patient, preventing retrograde flow

1. Data on file. © 2015 Quest Medical, Inc. Q2 and MultiPort are registered trademarks of Quest Medical, Inc. 2015-01


Distributed by:



997 Horan Drive | Fenton, MO 63026  
Customer Service: (800) 969-6331  
Email: [cs@progressivemedinc.com](mailto:cs@progressivemedinc.com)  
Fax: (314) 961-4535



# Q2® Split Septum Technology from QUEST Medical, Inc.



## ■ Controlling Infection Rates with Q2® Engineering

**“When needleless systems are used, a split septum valve may be preferred over some mechanical valves due to increased risk of infection with the mechanical valves.”** – CDC Guidelines for the Prevention of Intravascular Catheter-Related Infections; 2011:20, Department of Health & Human Services

Q2® Split Septum Needleless Connectors assist in meeting the OSHA Bloodborne Pathogen Standard by providing engineering controls which isolate or remove the hazard of bloodborne pathogen exposure from the workplace. These engineering controls coupled with work practice controls aid in eliminating or minimizing exposure which conform with paragraph (d)(2)(i) in OSHA's Methods of Compliance. 29 CFR 1910.1030, December 1991

Distributed by:



997 Horan Drive | Fenton, MO 63026  
Customer Service: (800) 969-6331  
Email: [cs@progressivemedinc.com](mailto:cs@progressivemedinc.com)  
Fax: (314) 961-4535



QUE9520revAUG16