



Vial2Bag® IV Admixture System

Universal, Needle-free Admixture System

The Vial2Bag IV admixture system enables safe and convenient reconstitution and transfer of a drug between a vial and an IV bag. The Vial2Bag system can assist in meeting the requirements of USP 797 for immediate use application (1 hour – Beyond Use Date) and can help optimize pharmacy clean room time, reduce drug and fluid waste and assist in standardizing admixture systems for nursing.

Vial2Bag IV Admixture System Features

- Compatible with both liquid and powder medications
- Conforms to any manufacturers' IV bags
- Swabable valve for easy connection to 13mm and 20mm vial adapters
- Needle free




Vial2Bag is ideal for

- Faster medication delivery time to nursing
- Waste reduction for medications and IV fluids
- Decreasing “choke points” in pharmacy
- Standardizing admixture systems on the nursing floor
- Assisting in meeting compliance with USP 797 immediate use

The Vial2Bag IV admixture system is 510(k) cleared by the United States Food and Drug Administration and carries the CE mark.

For more information visit www.westpharma.com

PROGRESSIVE MEDICAL

  	<p>Vial2Bag Item # 36098093</p> <p>13mm Orange Male Luer Vial Adapter Item # 36098094</p> <p>20mm Blue Male Luer Vial Adapter Item # 36098069</p>	<p>Quantity 50/Box</p> <p>150/Box</p> <p>150/Box</p>
--	--	---

Distributed by: Progressive Medical, Inc.
997 Horan Drive • Fenton, MO 63026

Ordering Information:
Phone: 314-961-5786 • Fax: 314-961-4535
cs@progressivemedinc.com