



**Distributed By**



Innovative Products. Uncommon Attention.

997 Horan Drive, Fenton, Missouri 63026-2401

P: 800-969-6331 or 314-961-5786

F: 314-961-4535

E: [cs@progressivemedinc.com](mailto:cs@progressivemedinc.com)

# GENICON<sup>®</sup>

## Nylon Specimen Retrieval






All-in-One Introducer Style Specimen Bags

## Laparoscopic Solutions: Specimen Retrieval Bag

The GeniStrong™ Bag is an easy-to-use retrieval system designed to temporarily contain specimen and facilitate its removal from the patient during laparoscopic surgery. Designed to minimize contamination of the abdominal cavity, the Specimen Retrieval Bag includes the following distinctive features:



### 550-000-008

2000mL  
14cm x 27cm  
15mm Introducer


### 550-000-005

750mL  
11cm x 22cm  
12mm Introducer

### 550-000-003

240mL  
7.3cm x 17cm  
10mm Introducer

## Burst Strength Comparison



\*Test Report Data On File



Bag-Only Specimen Bags



Length	Mouth Opening	Qty/Box	Volume	Part Number
3 in. / 7 cm	3.5 in. / 9 cm	20	50 mL	550-000-017
5 in. / 13 cm	3.5 in. / 9 cm	20	150 mL	550-000-009
6.5 in. / 17 cm	4 in. / 11 cm	20	240 mL	550-000-010
8 in. / 20 cm	8 in. / 20 cm	10	750 mL	550-000-011
10 in. / 23 cm	9 in. / 25 cm	10	1600 mL	550-000-012
13 in. / 34 cm	10 in. / 25 cm	10	3000 mL	550-000-013

# Clinically Preferred

Improved Patient Outcomes  
 Improved Patient Safety  
 Industry 1st Radiopaque Feature Allows Detection On X-rays  
 Improved Patient And Surgeon Satisfaction

# Cost Efficient

Available On Multiple GPO Contracts  
 Competitively Priced For The Health Care Provider  
 Reduced Procedural Costs  
 Can Be Used Multiple Times On A Single Patient

## Rip-stop Nylon

Rip-stop Nylon With Polymer Coating  
 For Leak-proof Performance  
 Improved Quality Of Care  
 Superior Strength Vs. Poly (Plastic) Bags

## Superior Strength

Superior Quality Provides Superior Performance  
 Does Not Burst Or Propagate Tears Like Poly (Plastic) Bags  
 Ideal For MIS Procedures  
 (Removing Large Specimen From Small Incisions)



### Recommended Anatomy Removal

Appendix, Ectopic Pregnancy, Gallbladder (Standard and Acute), Large Fibroid, Node Biopsy, Ovary, Prostate, Small Bowel Resection, Uterine Fibroids, VATS	550-000-003	550-000-010	550-000-000
Bowel Resection, Gastric Sleeve, Kidney, Lobe Resection (Liver/Lung), Spleen, Uterus, VATS	550-000-005	550-000-011	
Enlarged Spleen, Large Lobe Resection, Large Uterus, Malignant Kidney, Hemi-Colectomy, VATS	550-000-008	550-000-012 550-000-013	
Nanosurgery, Lymph Node Dissection, Oophorectomy, Pediatric Procedures	N/A	550-000-017 550-000-009	

### Features & Benefits

Feature	GENICON	Anchor	Applied	Medtronic	Ethicon
Rip-Stop Nylon Bag Material	✓	✓		✓	
Superior Strength	✓	✓			
Will Not Burst / Propagate Tears	✓	✓		✓	
Multi-Use (Single Patient)	✓	✓			
Available on U.S. GPO Contracts	✓		✓	✓	✓
Cost Efficient	✓		✓		
Multiple Sizes Available	✓	✓	✓	✓	
Gas Seal in Introducer	✓		✓		
Buttonhole Quick Bag Release Feature	✓				
Metal Pull/Push Rod (Sturdy Placement)	✓	✓	✓		✓
Radiopaque Feature Visible on X-Rays	✓				
Tapered Bag Shape for Easy Removal	✓	✓	✓	✓	✓
Translucent	✓	✓	✓	✓	✓